


⑤ lesser curvature of the middle body



<https://kotobukimedical.com/en/esd-training/G-master/setting/>



★ Training videos are available on the web page.

*For smooth setting, return the [X: Tube Left / Right] to the 0 position before setting.

*Be careful not to pinch your fingers when adjusting the [Z: Tube Up / Down].

Mark	Name	Location Number											
R	Mucosa Holding Arm Right	SET		1 (straight)		2 (curve)		S : 1 (gassing) ⇔ 2 (degassing)					
L	Mucosa Holding Arm Left	SET		1 (straight)		2 (curve)		3 (degassing)					
A	Vertical Rotation position	1	2	3	4	5	6	7	8	9	10	11	12
B	Horizontal Rotation Position	1	2	3	4	5	6	7	8	9	10	11	12
Y	Tube Front / Rear		+20	+10	0	-10	-20	-30	-40	-50			
X	Tube Left / Right	0	10	20	30	40	50	60	70	80	90	100	
			110	120	130	140	150	160	170	180	190	200	210
Z	Tube Up / Down	0	10	20	25	30	35	40	45	50	60	70	80
W	Tube Angle	-20	-10	0	5	10	15	20	25	30	35	40	40
M	Spatula Insertion Position	OFF	1	2	3	4	5						
	Back Rubber Strap	OFF	ON										